

Cloning 101

Adapted from TOPO TA Cloning Kit for Sequencing by Invitrogen (Life Technologies)

IMPORTANT

Before cloning: Make sure LB/AMP plates are made
No more than 8 samples to clone at one time

LB/AMP Plates:

- See recipe card or follow guidelines for grams/liter on container
 - Autoclave media and allow to cool before adding ampicillin
 - Add 0.1g ampicillin/liter
1. Begin with mixed sample
 2. Get genomic DNA (using zymo kit, MoBio kit, etc to get g DNA)
 3. Nanodrop
 4. Dilute to 50 ng/μl (with elution buffer or water)
 5. PCR 20 μl reaction (using NEB 2X MM)
 6. PCR clean-up kit
 7. Nanodrop
 8. Dilute to 10 ng/μl
 9. Time for cloning
 - a. Set LB/AMP plates out on bench to warm to room temperature
 - b. Get SOC media and competent cells from -80°C and x amount of tubes for the remaining samples (Be sure to take ice bucket to -80°C when retrieving competent cells)
 - i. Color lids purple and label with sample ID
 - ii. Set on ice
 - c. Set heat shock to 42°C and fill holes with water, according to number of samples or set water bath to 42°C
 - d. Get enough PCR tubes for number of samples (found under the PCR hood)
 - e. Get salt solution and water from -20°C
 - f. In each LABELED PCR tube: Add in the order shown

Reagent	Volume	Location
Water	Add to a total volume of 6 μl	-20°C
Salt Solution	1 μl	-20°C
Fresh PCR Product	0.5-4 μl	
TOPO Vector	0.5-1 μl	-20°C
Final Volume	6 μl	

Important! When working with vector, work efficiently in order to put it back in -20°C as quickly as possible

Keep all reagents on ice while working

- i. Before adding vector mix the water, salt solution, and DNA
 - ii. When adding vector, tap PCR tube carefully, do NOT mix vector
 - iii. Put samples at room temperature for 20 minutes for the ligation
- g. Separate competent cells for the number of samples
 - i. Example – the competent cells contain 50 μ l, so for 5 samples, separate 10 μ l in each purple colored tube. Don't split cells more than five ways
 - ii. Cut the tip off of the pipette tip using a razor blade so as not to damage the competent cells
- h. After incubation, place 2 μ l of ligation from PCR tube into LABELED purple colored tube, tap carefully (do NOT mix)
- i. Place on ice for 5-30 minutes
- j. Take ice bucket over to heat shock plate or water bath and heat shock for 30 seconds at 42°C
- k. Place immediately back on ice for 2-3 minutes
- l. Add 250 μ l SOC media to each tube, using a flame for sterility
- m. Tape tubes horizontally on a beaker and incubate while shaking for 1 hour at 37°C
- n. While incubating samples
 - i. Spread 50 μ l of x-gal onto LB/AMP plates then place in 37°C to warm up
- o. Put SOC media back in -80°C
- p. After incubation, get transformants from incubator, LB/AMP plates from 37°C, and hockey stick/ethanol/spinny device
 - i. Using a flame, put 60 - 75 μ l of transformants on plate and mix on LB/AMP plates with hockey stick device
- q. Place LB/AMP plates with transformants in 37°C overnight
- r. Put transformants tubes in beaker and store in refrigerator in case they're needed later