

BOIL and FREEZE DNA extraction protocol

1. Centrifuge the liquid culture, discard the supernatant
2. Froze the sample for 2 hours at -80°
3. Add $100\ \mu\text{l}$ of Buffer A and place the sample in a boiling water bath for 45 secs.
4. Add $30\ \mu\text{l}$ of lysozyme mix gentle by inversion. Incubate for 30 min at 37°C .
5. Add $10\ \mu\text{l}$ 20% SDS and $60\ \mu\text{l}$ proteinase K . Mix be gentle inversion. Incubate for 30 min. at $50\ ^{\circ}\text{C}$.
6. Add $100\ \mu\text{l}$ 20% SDS (approx 5% conc.) and $250\ \mu\text{l}$ phenol: CHCl_3 :IAA (24:24:1) ($245\ \mu\text{l}$ phenol and $255\ \mu\text{l}$ SEVAGE [CHCl_3 :IAA])
7. Spin at 12,000 g at $4\ ^{\circ}\text{C}$., pipette supernatant into a fresh 1.5 ml tube.
8. **IN HOOD FROM HERE:** Extract samples with 1 volume of phenol mix, spin at full speed for 2 min at 4°C . Remove aqueous layer and add 1 vol of CHCl_3 :IAA and mix (the DNA is with the aqueous layer. Spin at max speed for 2 min at 4°C . Remove the aqueous layer.
9. Precipitate nucleic acids by adding $24\ \mu\text{l}$ of 3M sodium acetate and $1000\ \mu\text{l}$ EtOH (search).
10. Centrifuge at max speed for 5 min. discard supernatant, Rinse pellet with $500\ \mu\text{l}$ of 70% ethanol. Air dry on kimwipe resuspend pellet in $15\ \mu\text{l}$ o elution Buffer of TE buffer.